


Configuration for Neuralynx Acquisition

The following scheme shows how is the Networks connections for the Neuralynx Acquisition and Saving Processes.



..... Private Network : 192.168.1. for MM-NLX Interface

----- Private Network : 192.168.0. for NAS connection

----- Clinic Network: Connection to the Clinic Network



MEDICAL EQUIPMENT

Private Network : 192.168.1.

for Micromed-Neuralynx Interface

This Network creates the private connection between MM Acquisition PC and Neuralynx PC. On it will pass only the data streaming of the acquisition from Neuralynx to Micromed Computer.

Each MM PC has the same configuration, similarly each NX PC has the same configuration: this is to avoid errors during the setting for acquisition.

MM PCs have the following IP address (on the second network card): 192.168.1.2; Subnet Mask: 255.255.255.0

NX PCs have the following IP address: 192.168.1.3; Subnet Mask: 255.255.255.0

Connection Procedure

1- Neuralynx side

Connect the Neuralynx on the Plug named NLX-NETCOM with a network cable. The Ethernet Interface must be the one with the label. If the light on the Card is blinking, the connection is ready. Open the Internet Explorer Browser and write :

\\192.168.1.2\System98\

If User and Password are required, insert

user: *aufnahme*

password: *eeg*

Once in the folder, go to the "Settings" subfolder and copy the file with the name "Micromed.cfg" to the Neuralynx Desktop. This file contains the settings for Cheetah program.

2- Micromed side

In SystemPLUS Software, select the following configuration on Extras→Systemeinstellungen→EEG:

Verstärker/Schnittstelle

- Digital Lynx Sx 256

- Network

Download the Connect and Montages files ready for Neuralynx Acquisition from the dedicated folder on the server, or prepare the connect files.

Once the two Computers are so set, the acquisition can start.

Note: in order to check the connections between the two PC's, open the command window on each PC and test with the command "ping -IP Address from the other PC-" and verify that the pinging run fine. If no, please check the network cable and the connection.

1. Strom - Computer an (zusätzliche Sicherung) + Verstärker
2. Netzwerk @ Verstärker an
3. Umbenennen des Channels (Micro-uAR) → RA
(R/L vom Anfang/Ende vertauschen) (~~Micro-RA~~)
4. Headstages anschließen (Referenz: Eingang/Rot
an Stirnelektrode Ausgang/Schwarz beliebig
(meist unter Schlüsselbein) Adapter
5. Cheetah Atlas starten (Dann umstellen lassen)
6. Head Boxen anschließen (2 Stück)

Micromed S.p.A.

Via Giotto 2, 31021 Mogliano Veneto, Treviso - Italy

Tel: +39-041-593 7000 - Fax: +39-041-593 7011

Partita IVA IT 03906850262

Reg. Impr. e C.F. TV 03906850262 - R.E.A. n. 307680



QUALITY SYSTEM CERTIFIED UNI EN ISO 9001:2008 UNI EN ISO 13485:2004

7. 2 Windows neuansetzen öffnen über Window → Add Window → Time
2 Macro 1 Obergliederung/ENG 8. rechte Maus klicken → zoom 8x

AR/RA
REC
RAH
RMH
RPHC
RAH
LA

Private Network : 192.168.0. ----- for NAS connection

This Network allows the Neuralynx Computers to connect to the NAS Systems in the Serverroom, through the Netgear Switch (it's allocated in Serverroom, Rack n.3).

NAS Configuration

Each DiskStation has a Volume, composed by 8 Disks in Synology Hybrid Raid (RAID 5) with a total capacity of 25,0 TB.

The screenshot shows the 'Volume' configuration page for 'Volume 1'. The 'Properties' section indicates it is a Synology Hybrid RAID (SHR) with a capacity of 25.04 TB and 1.07 GB used. The 'Disk Info' table lists 8 disks, all in 'Normal' status.

| Number | Disk Size | Status |
|--------|-----------|--------|
| Disk 1 | 3.64 TB | Normal |
| Disk 2 | 3.64 TB | Normal |
| Disk 3 | 3.64 TB | Normal |
| Disk 4 | 3.64 TB | Normal |
| Disk 5 | 3.64 TB | Normal |
| Disk 6 | 3.64 TB | Normal |
| Disk 7 | 3.64 TB | Normal |
| Disk 8 | 3.64 TB | Normal |

Each Volume has a Shared Folder named NLX_NAS_XX (where XX is the number of the NAS).


The screenshot shows the 'File Station' interface for the 'NLX_NAS_03' shared folder. The folder is empty, and the status bar at the bottom indicates 'Available: 25640 GB, Total: 25641'.

The first Network Interface is set:

Diskstation_1 : IP 192.168.0.21; Subnet Mask: 255.255.255.0

Diskstation_2 : IP 192.168.0.22; Subnet Mask: 255.255.255.0

Diskstation_3 : IP 192.168.0.23; Subnet Mask: 255.255.255.0




The second Network interface is connected to the Domain with this IP configuration:

Diskstation_1 : IP 10.47.11.111; Subnet Mask: 255.255.248.0

Diskstation_2 : IP 10.47.11.112; Subnet Mask: 255.255.248.0

Diskstation_3 : IP 10.47.11.113; Subnet Mask: 255.255.248.0



On the Domain you can see each NAS with the name:

Micromed S.p.A.
Via Giotto 2, 31021 Mogliano Veneto, Treviso - Italy
Tel: +39-041-593 7000 - Fax: +39-041-593 7011
Partita IVA IT 03906850262
Reg. Impr. e C.F. TV03906850262 - R.E.A. n. 307680



QUALITY SYSTEM CERTIFIED UNI EN ISO 9001:2008 UNI EN ISO 13485:2004

- DISKSTATION_1
- DISKSTATION_2
- DISKSTATION_3




Neuralynx Configuration

The second Network Card (the first is reserved for MM-NLX Interface) on the PC is set with a specific IP Address, in order to connect it automatically to the NAS System.




The IP Addresses are so set:

- Neuralynx 1 → IP 192.168.0.101; Subnet Mask: 255.255.255.0
- Neuralynx 2 → IP 192.168.0.102; Subnet Mask: 255.255.255.0
- Neuralynx 3 → IP 192.168.0.103; Subnet Mask: 255.255.255.0

Each PC has 3 different Network Drives already set:

-  nlx_nas_01 (\\192.168.0.21) (X:)
-  nlx_nas_02 (\\192.168.0.22) (Y:)
-  nlx_nas_03 (\\192.168.0.23) (Z:)

When available, the resources allow you to save data on them.

- ▶  nlx_nas_01 (\\192.168.0.21) (X:)
- ▶  nlx_nas_02 (\\192.168.0.22) (Y:)
- ▶  nlx_nas_03 (\\192.168.0.23) (Z:)




Connection Procedure

Connect the second Network Card of Neuralynx Computer to the Network Plug on the Wall. The name of the Plug is PRX-NLX-VTR.

If the light on the Network Card is blinking, the connection is ready.

Normally the user: *admin* and password: *admin* is required to access the resources.

Once the connection between Neuralynx Computer and NeteGear Switch is ready, directly on the computer Resources you can find the 3 NAS Resources mapped, with the name:

- ▶  nlx_nas_01 (\\192.168.0.21) (X:)
- ▶  nlx_nas_02 (\\192.168.0.22) (Y:)
- ▶  nlx_nas_03 (\\192.168.0.23) (Z:)


After this it's possible to save data on the remote folders.

Clinic Network : ----- to access to the data through a PC on the Network

You can access to the data saved on the NAS through the Clinic Network. Search on the Network Resources:

- \\DISKSTATION_1
- \\DISKSTATION_2
- \\DISKSTATION_3

User: *admin* and password: *admin* are required to access to the sheared folders.
Each Diskstation has one shared folder which can contain the data.



It's possible to map the network drive on the PC.

Micromed S.p.A.

Via Giotto 2, 31021 Mogliano Veneto, Treviso - Italy

Tel: +39-041-593 7000 - Fax: +39-041-593 7011

Partita IVA IT 03906850262


Reg. Impr. e C.F. TV03906850262 – R.E.A. n. 307680



QUALITY SYSTEM CERTIFIED UNI EN ISO 9001:2008 UNI EN ISO 13485:2004

Connection for Maintenance Processes


Connect the PC on the PRX-NLX-VTR in the room or directly on the NETGEAR Switch on the Serverroom.
Set the IP Interface on the Laptop with the IP: 192.168.0.80.
If you have installed Synology Assistant you can enter and manage the units directly with this tool.



If you don't have the program, insert directly on Internet Explorer:

- <http://192.168.0.21:5000/webman/index.cgi>
- <http://192.168.0.22:5000/webman/index.cgi>
- <http://192.168.0.23:5000/webman/index.cgi>


User: *admin* and password: *admin* are required in order to enter and manage the NAS Systems with Synology Assistant Software or Internet Interface.



Aquisition System Status → View → Acquisition Entity ^{Display Properties}
Macro: Headstage 001 / Subject Grand

Window 1, Window 2 (002)

(über die jeweiligen Time Windows wissen wir
welches Macro / Micro / Oberfläche Wires sind)

Micro: Headstage 005 Reference 1 oder 2 oder 3, 4
Oberfläche:  wie Macro
→ Input-Range auf $\pm 1000 \mu V$

Referenz: - hat am Stirn mit Verbindungs
Stecker

- Schwanze unter Schindelstein mit
Verbindungsstecker

- Grauer Endstück in RA 1

Neuen Ordner für Patienten mit
034 im Laufwerk E oder (F)
anlegen, dann Head-boxes anschauen
dann record

kommen auf Signal wegen Blutdruck-
messgerät / Lichtstrahler

- Settings ~~in~~ speichern in Config-file

- bei gutem Signal 58 record ^{mit dem Fenster} eingestellt werden
- Punkt → 17 no 11